

Hyperbaric Chamber Stand Type

Model: CH-HEL05 / CH-HEL06



Features:

1. 1.3 ATA (4psi) operating pressure.
2. 2 zipper seal for easier entrance.
3. Large transparent viewing windows to prevent claustrophobia.
4. Cotton chamber protection cover, to avoid dirty and easy to wash.
5. Emergency valve - To speed the depressurization in emergency.
6. User can observe the pressure while inside the chamber by internal pressure gauge
7. It can provide 93% oxygen from an oxygen headset/facial mask.
8. It is very easy for one people to operate without assistance.

Application:

1. Beauty Salon & Spa
2. Sports & Athletes & Gym
3. Clinic
4. Altitude Sickness
5. Home Use

Function of Hyperbaric Oxygen Therapy:

1. Reduce fatigue, improve sleeping quality and improve brain function.
2. Relieve work stress and support digestion.
3. Improve general wellness and increase energy level
4. Fast recovery from frequent flying or long flights
5. Fast recovery from sports injury,improve sports performance
6. Resolve alcohol concentration,regain consciousness after drunken
7. Improve cellular metabolism,circulation,inflation
8. Improve children's autism
9. Improve diabetic wounds
10. Improve cerebral palsy

www.hiperbarice.ro

HELVETIC IMPEX SRL

Romania | Brasov – 500058

Street Bucegi 1B

primo@sunjoy.ro

+40 755-098.486

Parameters:
1. Chamber CH-HEL05 / CH-HEL06

Pressure	1.3ATA(30kPa, 4psi)
Material	TPU
Size	130x160cm
Capacity	1-2 People


2. Chamber Bladder CH-HEL06


Pressure	1.3ATA(30kPa, 4psi)
Material	TPU
Size	180x175cm
Capacity	3-5 People


3. Air compressor

Voltage	110/220V
Oxygen Flow	72L/Minute
Power	480W
Weight	18kg


4. Oxygen Concentrator

Voltage	110/220V
Oxygen Flow	1-5L/Minute
Power	480W
Weight	30kg


5. Air Cooler

Voltage	110/220V
Capability	150L/Minute
Pressure	0.8mPa
Power	480W
Weight	30kg


10. Improve cerebral palsy

www.hiperbarice.ro

HELVETIC IMPEX SRL

Romania | Brasov – 500058

Street Bucegi 1B

primo@sunjoy.ro

+40 755-098.486